**Manufacturer:**  
Gebr. Brasseler GmbH & Co. KG  
Trophagener Weg 25  
32657 Lemgo · Germany

**Distributor/Importer for USA and Canada:**  
Komet USA LLC  
3042 Southcross Blvd, Suite 101  
Rock Hill, SC 29730, USA

Phone: 888-566-3887  
Fax: 800-223-7485  
info@kometusa.com  
www.kometusa.com

**CAUTION: USA and Canada:**  
devices to sale by or on the order of a licensed dentist.

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**Instructions for use**  
Please read carefully and retain for future reference.

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**General Instructions for use and safety recommendations for the application of rotary and oscillating dental instruments**

Rotary and oscillating instruments for use in the dental sector are to be used only by doctors or other experts who are experienced in the safe handling of these instruments.

**Area of application**  
These general instructions for use and safety recommendations apply to all products and are to be generally observed. If separate instructions for use are provided in the packaging, these take precedence over these general instructions for use.

**Storage**  
Protect packed instruments from UV rays and high temperature. Store all instruments in a clean and dry environment. Do not store in the same room as solvents or chemicals. Instruments without direct identification marking have to be stored with their original packaging.

- 1. Proper use**
- Make sure that only technically and hygienically perfect and cleaned turbines, handpieces and contra-angles are used.
  - Chuck the instruments as deeply as possible.
  - The instrument must be rotating before contact is made with the tooth.
  - Medical devices which are identified for single use (marked on the packaging) are not intended for reprocessing. Recommendations for safe reprocessing are therefore not supplied.
  - Avoid jamming and using the instrument as a lever as this leads to an increased risk of fracture.
  - Wear safety glasses as required.
  - Avoid unprotected contact with the instruments (use protective gloves).
  - Thermal damage caused by rotary instruments has to be avoided in any case (work at recommended speed and use sufficient water cooling).
  - Coarse and super coarse grit diamond instruments (including S-Diamonds®, Series 2000 and Turbo-Diamonds) may lead to increased thermal stress. Therefore, when using such products, use sufficient water cooling (at least 50 ml/min) and work at minimal contact pressure. To achieve an optimal surface quality, subsequent finishing is necessary.
  - Preferably use instruments with rounded edges as the preparation of sharp-edged undercuts may lead to an increased risk of a damaging notch effect. Improper use leads to increased risk and inferior results. Therefore, stick to the application and speed recommendations indicated on the labels and in our instructions for use.

- 2. Recommended speeds**
- The general rule is:
- The larger the working part, the lower the speed
  - Maximum speed  $\varnothing_{max}$  300 000 rpm means: Suited for micromotor handpieces and turbines with stable ball bearings. Not recommended for old turbines with air bearing.
  - Maximum speed  $\varnothing_{max}$  30 000 - 160 000 rpm means: Suited for micromotor handpieces or lab handpieces up to the speed indicated. Not recommended for turbines. Not observing the maximum permissible speed leads to an increased safety risk.

- 3. Contact pressure**
- Excessive contact pressure (>2N) has to be avoided.
- with cutting instruments, this can lead to damage to the working part and to chipping of the blades as well as an excessive generation of heat.
  - Increased contact pressure may lead to stripping of the grit on abrasive instruments or to clogging of the instruments and increased heat generation.
  - Increased contact pressure may also lead to thermal damage to the pulp or, in case of damaged blades, to rough surfaces. In the extreme cases, instrument breakage may even occur.

- 4. Cooling**
- To avoid undesirable heat generation during preparation, make sure to provide sufficient cooling by means of air/water spray (at least 50 ml/min.).
  - Additional external cooling is required when using FG instruments with a total length of more than 22 mm or a head diameter exceeding 2 mm.
  - Insufficient water cooling can result in irreversible damage to the tooth and the surrounding tissue.

- 5. Guideline on the number of times rotary instruments can be used**
- The below values are guidelines. The service life of the instruments may differ from these values as this depends on the application and/or the material treated.
- In certain cases, the instruments can be used more often, provided that there are no visible signs of wear.

Single Use Instrument:	- 1 x
Instruments made of unprotected tool steel:	- 4 x
Instruments made of stainless steel:	- 4 x
Instruments with internal cooling:	- 4 x
Tungsten carbide instruments:	- 15 x
Diamond instruments:	- 25 x
Polishers:	- 10 x
Ceramic abrasives:	- 10 x
Endodontic instruments - wide canals:	- 8 x (max.)
Average canals:	- 4 x (max.)
Narrow canals:	just use 1 x
For hygienic reasons, polishers with lamellae and dental brushes may only be used	- 1 x

**Disposal**  
To prevent contamination, discard instruments in tight, puncture resistant containers.

- 6. Elimination of worn instruments**
- Komet® Single Use Instruments are pre-sterilized for single-use and should not be reprocessed or reused.
  - Damaged and deformed blades cause vibrations and lead to poor preparation margins and rough surfaces.
  - Blank spots on the surface of diamond instruments are an indication of abrasive grit wear and reduced cutting efficiency. These deficiencies lead to excessive temperature and finally pulp damage. Therefore, bent or "out of true" instruments must be eliminated immediately.
  - Very important: Blunt and damaged instruments may lead to the dentist applying higher contact pressure which may result in an increased operating temperature and thermal damage to the pulp. Damaged instruments therefore have to be discarded immediately.

- 7. Sterilization procedure**
- Instruments that are supplied non sterile have to be prepared prior to first use.
  - Instruments that are supplied sterile do not require sterilization prior to first use.
  - Prior to the first use on the patient and immediately after each use, all rotary instruments have to be (re)processed.
  - Wear puncture resistant utility gloves when handling contaminated instruments.
  - Pre-soak instruments immediately after use to loosen debris.
  - For cleaning we recommend to use automated cleaning equipment (e.g. ultrasonic cleaner or washer-disinfector).
  - For manual cleaning use an ultrasonic cleaner with an enzymatic detergent (e.g. Cidezyme). For manual disinfecting use an ultrasonic cleaner with a disinfecting detergent (e.g. Cidex OPA).
  - For mechanical reprocessing a washer-disinfector with an alkaline detergent (e.g. Neodisher Mediclean forte)
  - Burs should be separated from each other in a bur block to prevent damage during immersion.
  - Brush away remaining debris with a metal cleaning brush (9791) or - especially in case of cleaning ceramic instruments - with a softer nylon cleaning brush (9873) and rinse instruments under running water.

- Dry burs (e.g. by airblasting) to avoid corrosion.
- Sterilization of instruments is carried out in the autoclave using the known methods. Minimum hold times: for 4 minutes at 132°C and 20 minutes drying. Times are hold times, running times are longer and may vary from unit to unit. The instructions provided by the device manufacturer have to be observed.
- After sterilization, check instruments for defects on the surface (corrosion). Do not reuse corroded instruments.
- The operator of medical products is responsible for seeing that proper treatment is carried out by qualified staff using the appropriate materials and suited equipment.
- For polishers, brushes and instruments with internal cooling please refer to the following specific information.

- 8. Specific instructions for individual instrument types**
- Tungsten carbide**
- When trimming dry plaster a suction device must be used.
  - Avoid any contact with H<sub>2</sub>O<sub>2</sub> (hydrogen peroxide). The carbide working parts would be attacked and damaged reducing the instrument's working life.

- Steel**
- Tool steel instruments may corrode when sterilized in the autoclave (see information in the catalog table).
  - When using the separating strip (9816) please avoid contact with the gingiva as there is a risk of injury due to sharp blades.

- Diamond**
- Use a disc guard for rotating diamond discs when working intraorally.
  - Use diamond discs 8964 and 983 only in standard clockwise rotation. Counter-clockwise rotation may lead to jamming and breakage.
  - When using finishing disc 952 avoid axial deflection over 45° and radial deformation. This may lead to breakage. In order to guarantee non-traumatic treatment, ensure direct vision and avoid contact with soft tissue. An integrated sliding clutch stops the disc in case it gets jammed. Afterwards the finishing disc can no longer be used.
  - When using the diamond strips (numbers start with WS or DS) please avoid contact with the gingiva as there is a risk of injury. Please also avoid extreme bending as this might lead to fracture of the strip.
  - Coarse and super coarse grit diamond instruments (including S-Diamonds®, Series 2000 and Turbo-Diamonds) may lead to increased thermal stress. Therefore, when using such products, use sufficient water cooling (at least 50 ml/min) and work at minimal contact pressure. To achieve an optimal surface roughness, subsequent finishing is necessary.

- Ceramics**
- When using ceramic instruments for cutting, care should be taken that the instrument does not get jammed; this increases risk of fracture.
  - Do not use the instrument as a lever.
  - Only use brushes with metal-free bristles for pre-cleaning (9873).

- Instruments featuring internal cooling**
- Special care has to be taken when cleaning the long and narrow bores and blind holes of these instruments.
  - Information with regard to cleaning: Clean the workplace and surfaces from all contaminations with a single-use cloth or tissue. Preparation of instruments with internal cooling has to take place within one hour after they have been used at the latest. It is not necessary to disassemble the instruments with internal cooling. Cleaning is carried out by means of a mandrel (9793), a cleaning brush (9791), a 10 ml syringe and running water.
  - Manual cleaning: Penetrate the internal cooling channel with the mandrel until it is clean without any further contamination at all. Brush the surface of the instrument under running water until all contamination has disappeared. Subsequently, the instruments with internal cooling have to be rinsed thoroughly by means of the 10 ml syringe filled with demineralised water.

- Bone cutters (tungsten carbide, stainless steel, diamond, ceramics)**
- Bone cutters are suitable for a wide range of uses in dental alveolar surgery, depending on the shape of their working parts, for example: for the restoration of edentulous ridges when removing sharp bone edges, for the extraction of bones for example from the chin or from the retromolar zone for autogenous bone transplantation, for osteotomy when exposing impacted teeth and for the treatment of root apexes as surgical measure for the conservation of teeth. Recommended speeds are mentioned on the label of the package of an instrument.
- Special care has to be taken when using the H254E/LE or H162SL, minimally invasive combination instruments for conservative preparation of bone tissue and hard tooth substance.
  - The H254E/LE and H162SL (FG shank) are designed to be used in the micromotor (red contra-angle) and not in the turbine → risk of accident.

- Trepan burs (stainless steel)**
- Trepan burs are indicated for removal of endosseous implants (227A, for implants of 18 mm length maximum), or for the extraction of bone transplants from suited donating areas in the oral cavity (227B, TRC, TRT).
- In order to eliminate any danger to the adjacent structures such as tooth root, tooth root, maxillary antrum or nerve canal the area surrounding the extraction site has to be carefully examined.
  - Avoid leverage of the trepan bur. Otherwise there is the risk that the trepan bur may slip off of the intended extraction site and cause injury.
  - Observe the optimum speeds on the labels, drill with constant exterior cooling with light pressure.
  - A sterile procedure is important to avoid infections.

- Root canal reamers (stainless steel)**
- Gates Glidden Bur Type "G", G180, G180A
- Green contra-angle,  $\varnothing_{max}$  450 - 800 rpm. For preparation of the coronal portion of the root canal, before or after the use of files or K-burs.
- Pulp Bur round 191
- Green contra-angle,  $\varnothing_{max}$  450 - 800 rpm. For root canal preparation.
- Reamer "P" 183L
- Green contra-angle,  $\varnothing_{max}$  800 - 1200 rpm. For preparation of the root canal entrance and enlarging the coronal portion, for the removal of gutta percha.

- Anodized Aluminum Bur Block**
- Prior to sterilization, rinse bur block under running water and make sure, that especially the holes are properly dried (e.g. by airblasting).
  - For the sterilization of instruments it is necessary to seal the bur block in sterilization pouch and sterilize it in the autoclave.
  - Multiple sterilizations of the bur block may lead to minor color deviations. Such a change in color does not have any effect on the quality or the performance of the bur block.
  - Due to its material the bur block is neither suitable for the sterilization in cold sterilization solutions nor for the washer/disinfector.

- Root posts**
- Root posts made of fiber-reinforced composite (ER DentinPost®, DentinPost®X), made of zirconium oxide ceramics (ER C-Post™) and made of polymethylmethacrylate (ER CAST) cannot be sterilized. Therefore, they must be cleaned and disinfected with medical alcohol. Root posts are intended for single use only.

- 9. Range of application of abrasives**
- High-grade corundum**
- Brown abrasives are intended for rapid grinding of metal alloys - Coarse grinding
  - Pink abrasives are for universal grinding of metal alloys - Medium grinding
  - White abrasives retain their original shape during grinding - Fine finishing

**Silicon carbide**  
- Green abrasives are intended for grinding ceramics  
**Athermon Discs**  
- Abrasive Discs intended for grinding ceramics  
- Recommended speed: ⌀<sub>int</sub> 15 000 rpm (⌀<sub>max</sub> 20 000 rpm!)  
- Apply low contact pressure (max. 2N), avoid jamming and leverage  
- Do not allow products to come in contact with humidity: Store and grind dry, do not use liquid coolant  
Attention: Overaged Athermon Discs may break spontaneously!  
- Eye protection and working with a suction system are recommended.  
- Dress Discs with suitable tools, e.g. Komet® Dressing Diamond 150.16.000, using 15 000 rpm with low contact pressure

#### 10. Maintenance of diamonds with sintered bond (DSB)

For cleaning and sharpening, DSB abrasives have to be dressed with the white cleaning stone 9750 from time to time:  
- Wet cleaning stone to avoid dust and trim the DSB abrasive to be cleaned/sharpened at ⌀<sub>int</sub> 15 000 rpm.  
- It is recommended to subsequently clean the trimmed DSB abrasive in the ultrasonic bath for 1- 2 minutes.  
- To exclude any damage to the instrument shank, avoid contact of the cleaning stone with the shank during the grinding procedure.

Attention! Bond abrasion can cause discoloration during reduction of ceramic restorations with sintered diamond tools. To remove any discoloration we recommend a brief final touch-up using green stones or sandblasting the surface.

#### 11. Polishers/Brushes

- Apply low contact pressure in order to minimize heat generation.  
- Apply a speed of ⌀<sub>int</sub> 5 000 - 6 000 rpm.  
- Polishing should always be carried out in circular motion.  
- In order to achieve a high shine polish, for multiple step polishing systems all polishers are to be used in the indicated sequence.  
- Use breathing mask (mouth and nose) as well as a suction device in the laboratory.  
- Eye protection is recommended.

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1. Désinfectant and cleaning  
Due to their material properties, brushes and polishers have to be cleaned differently from rotary instruments. Use disinfecting and cleaning agents that are suited for polishers. Mix the solution by observing the indicated concentration.  
2. Sterilization  
Autoclave only.

#### 12. Safety and liability

See also additional information in the Komet USA LLC catalog and our website at [www.komet-usa.com](http://www.komet-usa.com). The user is responsible to check the products prior to use to determine whether they are suited for the intended purpose. The instruments should only be used according to their intended application.

In case of contributory negligence by the user, the manufacturer partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

#### 13. Disclaimer of Warranties

ALL WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXPRESSLY DISCLAIMED. ALL INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO, LOST PROFITS, INJURY TO OTHER PROPERTY, LOSS OF USE OR OTHER COMMERCIAL LOSSES, ARE EXPRESSLY DISCLAIMED. If, not withstanding the foregoing, consequential and incidental damages cannot be excluded due to operation of law, such damages are expressly limited in amount to the purchase price of any defective products.

Store out of children's reach.  
For dental use only.

It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. The user is responsible for the application of the instruments. The manufacturer declines any liability for possible damage caused in connection with the use of this instrument.

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© Gebr. Brasseler GmbH & Co. KG